# ArthroZheal

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# ArthroZheal

## Introduction

ArthroZheal is an autologous platelet-rich fibrin product produced by the Vivostat® system.

The Vivostat® system comprises: A Processor Unit, a disposable Preparation Kit for preparation of ArthroZheal, an Applicator Unit with a footswitch and a disposable Application Kit used to apply ArthroZheal during surgery.

Optionally, a Co-delivery Kit enables applying an additional liquid substance (antibiotics, stem cells, chondrocytes, etc.) with ArthroZheal via a separate lumen in the tubing line. Up to 5.3 ml of substance may be applied together with ArthroZheal. The choice of an additional liquid substance and its influence on the patient is the responsibility of the treating physician.

This Preparation Kit/Application Kit is part of a system and these instructions should be read in conjunction with the Vivostat® user Manuals and ArthroZheal Quick Reference Cards.

**ArthroZheal Endoscopic Handle:** The ArthroZheal Endoscopic Handle is a reusable device, used to insert and direct the endoscopic catheter.

## Intended Use

The Vivostat® system is a medical device used for the preparation and application of ArthroZheal from whole blood.

ArthroZheal has different scopes of application: Tissue regeneration, haemostasis, tissue sealing and tissue gluing.

## PREPARATION KIT – INSTRUCTIONS FOR USE Preparation Kit contents

A sterile disposable single-use Kit, which contains all the components needed for the preparation of ArthroZheal.

● **Preparation Unit** – The Preparation Unit is a device that extracts platelets and fibrin from the blood and thereby develops ArthroZheal. A tubing line is attached to the Preparation Unit used for transferring the blood into the Preparation Unit.

● **Citrate/TA (ArthroZheal)** – A vial containing a sterile solution of citrate (anticoagulant) and 300 mg tranexamic acid (antifibrinolytic).

● **Spike** – The spike is connected to the tubing line and used to pierce the Citrate Vial.

● **16 G Needle** – The 16 G needle is connected to the tubing line of the Preparation Unit and is used to draw blood from the patient.

● **Sugi Swab** – The two non-sterile sugi swabs are used to clean the valve after the tubing line has been removed from the Preparation Unit.

● **pH4** – The pH4 syringe is loaded into the Preparation Unit prior to processing. **After processing, this syringe contains the platelet-rich fibrin solution and is referred to as the ArthroZheal syringe.**

● **Replacement Cap** – The Replacement Cap is used to protect the luer tip of the ArthroZheal syringe.

● **Disinfection Swab** – Used for disinfection of the rubber membrane on the Citrate/TA.

● **Removable Batch Tracking labels (4 pcs.)** – Used for documentation in the patient journal.

## Filling the Preparation Unit with Blood

ArthroZheal can only be prepared from patient-derived donation of blood. We recommend drawing blood directly from the patient before any procedure is initiated.

Please notice! When the Preparation Unit contains the blood it must be handled carefully. Avoid shaking movements from side to side. The Preparation Unit containing the blood must be processed within one hour after collecting the blood.

After preparation, immediately remove the Preparation Unit from the Processor Unit, turn it upside down and remove the ArthroZheal syringe from the Preparation Unit.

After preparation, store the ArthroZheal syringe at room temperature. The ArthroZheal syringe should be used within eight hours of preparation and be carefully turned upside down 10 times to mix the platelets with the fibrin right before use.

## APPLICATION KIT – INSTRUCTIONS FOR USE Application Kit contents

The Application Kit is available as a Spraypen or an Endoscopic Applicator

A sterile disposable single-use kit, which contains all the components needed to apply ArthroZheal.

● **Spraypen/Endoscopic Applicator** – The application device is used to apply ArthroZheal.

● **pH 10** – The pH 10 solution activates the polymerisation of the ArthroZheal upon spraying.

● **pH Paper** – The pH paper is used to test adequate mixing of pH10 and ArthroZheal.

● **1 ml Syringe** – The 1 ml syringe is used to deliver pH 10 buffer via the application device.  
Caution: Avoid air bubbles when filling with pH 10

● **21 G Needle** – The 21 G needle is used to draw pH10 buffer into the 1 ml syringe.

● **Batch Tracking labels (4 pcs.)** – Used for documentation in the patient journal.

## Contraindications

Known hypersensitivity to batroxobin and tranexamic acid.

ArthroZheal is not to be used in neurosurgery due to risk of hyperexcitability induced by tranexamic acid.

## Cautions and Warnings

**Patient-derived ArthroZheal is for autologous use only – check patient details to ensure patient donor and recipient are the same.**

Spray is achieved using pressurised air. Caution should be used where there is a risk of air embolism. We recommend using the spray mode "Jet No Air" in endoscopic applications.

Not for intravascular use.

Dispose of properly after use.

Universal precautions for blood handling should always be used in the operation of this product. Dispose of all accessories that may contain biohazard materials, such as body fluids, using universal blood handling precautions.

### **Interactions**

Drugs that reduce fibrinogen levels such as anicrod may affect the preparation of the ArthroZheal. **Administration of cefuroxime just prior to drawing blood may affect the coagulation abilities of ArthroZheal.** The effect of acute thrombolytics such as streptokinase, urokinase and t-PA on the preparation and fibrinolytic stability of the ArthroZheal has not been evaluated and cannot be excluded.

Studies have been conducted to demonstrate that heparin and warfarin therapy do not affect the preparation of ArthroZheal using the Vivostat® system.

Platelet aggregation and release of platelet-derived growth factors may be affected for up to 7 days after ingestion of analgesics containing acetylsalicylic acid. Other anti-platelets drugs (analgesics) may also have the same effect. The consequence of this for performance of ArthroZheal is not known currently.

Thrombocytopenia in the patient may cause lower contents of platelets than expected.

No other interactions are currently known.

### **Packaging, Storage and Shelf Life**

Store between 4°C and 25°C.

Do not use after the expiry date on package labels.

Do not use if the packaging of any components appears open or damaged.

### **Product and Patient Tracking**

A patient ID label should be placed on the ArthroZheal syringe and not on the Preparation Unit, as this will interfere with the process.

The label placed on the ArthroZheal syringe should be checked prior to use for patient identification purposes.

### **Application Technique**

- The application of ArthroZheal can only be activated with the Vivostat® footswitch.
- The spray mode can be adjusted using the button on the Applicator Unit.
- In case of blockage of the application device, cut off up to 8 mm of the tip of the application device with a sharp scalpel - do not use scissors.
- If you change to a new application device always remember to refill the pH10 syringe.

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Vivostat A/S