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VIVOSTAT® SYSTEM

Introduction

This Fibrin-Preparation Kit / Co-Delivery Kit is part of a system and these instructions should be read in conjunction with the Vivostat® User Manuals and Quick Reference Cards.

The Vivostat® system comprises a Processor Unit, a disposable Fibrin-Preparation Kit for preparation of fibrin, an Applicator Unit with a footswitch and a disposable Co-Delivery Kit used to apply the sealant and the desired substance (drug, stem cells etc.) during surgery.

Intended use

The Vivostat® system is a medical device used for the preparation and application of a fibrin sealant from a whole blood or plasma source. The Co-Delivery Kit has the ability to co-apply a fluid (substance) in a separate tube. Choice of substance and its influence on the patient is the surgeon's responsibility.

Fibrin sealants have been shown to have a range of uses including haemostasis, tissue sealing (including fluid and air) and tissue gluing.

FIBRIN PREPARATION KIT – INSTRUCTIONS FOR USE

Vivostat® Fibrin Preparation Kit Contents

A sterile disposable single use kit, which contains all the components needed to prepare fibrin:

- **Fibrin Preparation Unit** - The Fibrin-Preparation Unit is a device where fibrinogen in the blood/plasma is converted to fibrin. A Tubing Assembly is attached to the Fibrin-Preparation Unit used to fill the Fibrin-Preparation Unit.
- **Citrate/TA** - a vial containing a sterile solution of citrate (anticoagulant) and 100 mg tranexamic acid (antifibrinolytic).
- **Codan Spike** - The dispensing spike is used to access the citrate vial and connected to the Fibrin-Preparation Unit.
- **Sugi swab** - The two non sterile Sugi swaps are used to clean the valve after the tubing assembly has been removed from the Fibrin-Preparation Unit.
- **pH4** - the pH4 syringe is loaded into the Fibrin-Preparation Unit prior to processing. After processing, this syringe contains the fibrin solution and is referred to as the fibrin syringe.
- **Replacement Cap** - The Replacement Cap is used to protect the luer tip of the fibrin syringe.
- **Disinfection swab** - Used for disinfection of the rubber membrane on the Citrate/TA.
- **Batch Tracking labels (4 pcs.)** - used for documentation in the patient journal.

Filling the Fibrin Preparation Unit with Blood/Plasma

You can fill the Fibrin Preparation Unit up to 24 hours before you place it in the Processor Unit. If the blood is drawn more than four hours before use in the Processor

Unit, the Preparation Unit should be kept in the refrigerator at 5 degrees C.

The Vivostat® sealant may be prepared either from a patient-derived peri-operative donation, from pre-donated blood/plasma or blood from the heart-lung machine.

Pre-donated Blood/Plasma (anticoagulated bag) or blood from heart-lung machine

We recommend adding tranexamic acid as an antifibrinolytic agent by adding an approved sterile injectable solution of tranexamic acid (e.g cyklokapron® injection 100mg/mL). Alternatively the citrate vial supplied with the kit can be used, but this will reduce the amount of fibrin (approximately 20%).

1. Connect a 3-way valve to the Luer Connector. Add 100mg tranexamic acid sterile solution to the Fibrin Preparation Unit. Now remove the 3-way valve.
2. Connect the luer connector on the inlet line to the blood bag or the heart-lung machine and allow blood to flow into the Fibrin-Preparation Unit. Warning: Heavy blood flow from heart-lung machine might cause bubbles.

After preparation, immediately remove the Fibrin Preparation Unit from the Processor Unit, turn it upside down and remove the fibrin syringe from the Fibrin Preparation Unit.

The CO-DELIVERY KIT - INSTRUCTIONS FOR USE

This Co-Delivery Kit is part of a system and these instructions should be read in conjunction with the Vivostat® User Manuals and Quick Reference Cards.

The Co-Delivery Kit only works together with the Applicator Unit APL 404.

The co-application of fibrin and additional substance can only be activated with the Vivostat® footswitch as there is no pen button on the Co-Delivery Spraypen.

Vivostat® Co-Delivery Kit Content

A disposable single-use kit, which contains all the components needed to apply the sealant.

- **Co-Delivery Spraypen/Endoscopic Catheter** - Used to mix and apply the fibrin with the additional substance.
- **pH Paper** - The pH paper is used to test adequate mixing of the fibrin.
- **pH 10** - The pH 10 solution activates fibrin polymerisation when spraying.
- **21 G Needle** - The 21 G needle is used to draw the pH10 Buffer into the 1mL syringe
- **1 ml Syringe** - The 1 ml syringe is used to deliver pH 10 buffer via the application device.
- **Co-Delivery Syringe** - For delivering a fibrin compatible substance via the Co-delivery Spraypen. Caution: Avoid air bubbles in both syringes when filling with pH-10 and the fibrin compatible substance. N.B. You can not fill more substance in the Co-Delivery syringe than is found in the fibrin syringe. See illustration on the last page.

• **Replacement Cap** - The Replacement Cap is placed on the tip of the Co-Delivery syringe if it is not used immediately after filling.

• **Batch Tracking labels (4 pcs.)** - used for documentation in the patient journal.

Contraindications

Known hypersensitivity to batroxobin and tranexamic acid.

Known hypersensitivity to the chosen additional substance.

Cautions and Warnings

It is recommended that, whenever possible, Vivostat® Sealant should be prepared from patient-derived blood or plasma.

Patient-derived sealant is for autologous use only – check patient details to ensure patient donor and recipient is the same! Discard any remaining sealant at the end of the procedure.

Note: When using blood or plasma derived from a source other than the patient, the risk of transmission of blood-borne infectious agents cannot be entirely excluded.

Spray is achieved using pressurised air. Caution should be used where there is a risk of air embolism.

Direct topical application to the central nervous system of sufficiently high amounts of tranexamic acid (including application in fibrin sealants) may induce hyperexcitability, which may lead to convulsions. However, the concentration of tranexamic acid in Vivostat® (about 0.5mg/mL) has been shown to be between 6 and 12 times lower than that required to produce any noticeable effects in a rat study, where it was applied directly to the exposed spinal cord. Furthermore, clinical studies using Vivostat® in neurosurgery have identified no adverse effects related to the use of Vivostat®. However, the means by which hyperexcitability is induced by tranexamic acid and the importance of other factors when applied in fibrin sealants (e.g. individual variation, threshold in concentration/amount applied, formulation and release characteristics, mode and site of application) are poorly understood. Therefore, although indications are that the risk of inducing convulsions with Vivostat® is small, the possibility of inducing convulsions in susceptible individuals cannot be ruled out. Therefore, where direct contact of the Vivostat® solution is made with the central nervous system, such as in spinal or neurosurgery, the operating staff should be aware of the possibility of convulsions and be prepared to take appropriate corrective action if they should occur.

Choice of additional substance and its influence on the fibrin is the surgeon's responsibility.

Not for intravascular use.

Dispose properly after use.

Universal precautions for blood handling should always be used in the operation of this product. Dispose of all accessories that may contain biohazardous materials, such as body fluids, using universal blood handling precautions.

Interactions

Drugs that reduce fibrinogen levels such as batroxobin and anicrod may affect the preparation of the sealant. **Administration of cefuroxime just prior to drawing blood may affect the coagulation abilities of the fibrin.** The effect of acute thrombolytics such as streptokinase, urokinase and t-PA on the preparation and fibrinolytic stability of the sealant has not been evaluated and cannot be excluded. Studies have been conducted to demonstrate that heparin, warfarin and aspirin therapy do not affect the preparation of fibrin using the Vivostat® System.

No other interactions are currently known.

Packaging, Storage and Shelf Life

Store between 4°C and 25°C.

Do not use after the expiry date on package labels

Do not use if the packaging of any components appears to be open or damaged.

After preparation, store the fibrin syringe at room temperature. The fibrin syringe should be used within eight hours of preparation.

Product and Patient Tracking

Labels should be placed on the fibrin syringe and not on the Fibrin Preparation Unit, as this will interfere with the process.

Labels placed on the fibrin syringe should be checked prior to use for patient identification purposes.

Application Technique

- Dab the area of application until it is dry.
- The spray mode on the APL 404 Applicator can be adjusted on the User Button. We recommend to use the spraymode „low“ for Co-Delivery applications.
- In case of blockage of the Spraypen, cut off approximately 8mm of the tip of the spraypen with a sharp scalpel.
- If you change to a new Spraypen always remember to refill the pH-10 syringe.

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Vivostat A/S