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THE VIVOSTAT® SYSTEM

Introduction

This Fibrin Preparation Kit / Spraypen Kit / Endoscopic Kit is part of a system and these instructions should be read in conjunction with the Vivostat® User Manuals and Quick Reference Cards.

The Vivostat® system comprises a Processor Unit, a disposable Fibrin Preparation Kit for preparation of fibrin, an Applicator Unit and a disposable Application Kit used to apply sealant during surgery.

Endoscopic Applicator Handle: The Endoscopic Applicator Handle is a reusable device, used to insert and manipulate the Endoscopic Applicators.

Intended use

The Vivostat® system is a medical device used for the preparation and application of a fibrin sealant from whole blood or a plasma source.

Fibrin sealants have been shown to have a range of uses including haemostasis, tissue sealing (including fluid and air) and tissue gluing.

FIBRIN PREPARATION KIT - INSTRUCTIONS FOR USE

Vivostat® Fibrin Preparation Kit contents

A sterile disposable single use kit, which contains all the components needed to prepare fibrin.

- **Fibrin Preparation Unit** - The Fibrin Preparation Unit is a device where fibrinogen in the blood/plasma is converted to fibrin. A tubing assembly is attached to the Fibrin Preparation Unit.
- **Citrate/TA** - a vial containing a sterile solution of citrate (anticoagulant) and 100 mg tranexamic acid (antifibrinolytic).
- **Spike** - The dispensing spike is used to access the citrate vial and connect to the Fibrin Preparation Unit.
- **Sugi swab** - The two non sterile Sugi swaps are used to clean the valve after the tubing assembly has been removed from the Fibrin Preparation Unit.
- **pH4** - the pH4 syringe is loaded into the Fibrin Preparation Unit prior to processing. **After processing, this syringe contains the fibrin solution and is referred to as the fibrin syringe.**
- **Replacement Cap** - The Replacement Cap is used to protect the luer tip of the fibrin syringe.
- **Disinfection swab** - Used for disinfection of the rubber membrane on the Citrate/TA.
- **Removable Batch Tracking labels (4 pcs.)** - Used for documentation in the patient journal.

Filling the Fibrin Preparation Unit with Blood/Plasma

You can fill the Fibrin Preparation Unit up to 24 hours

before you place it in the Processor Unit. If the blood is drawn more than four hours before use in the Processor Unit, the Preparation Unit should be kept in the refrigerator at 5 degrees C.

The Vivostat® sealant may be prepared either from a patient-derived peri-operative donation, from pre-donated blood/plasma or blood from the heart-lung machine.

Pre-donated blood/plasma (anticoagulated bag) or blood from heart-lung machine

We recommend adding tranexamic acid as an antifibrinolytic by addition of an approved sterile injectable solution of tranexamic acid (e.g. Cyclokapron® injection 100mg/ml). Alternatively the citrate vial supplied with the kit can be used, but this will reduce the amount of fibrin (approximately 20%).

1. Connect a 3-way valve to the luer connector. Add 100 mg tranexamic acid sterile solution to the Fibrin Preparation Unit. Now remove the 3-way valve.
2. Connect the luer connector on the inlet line to the blood bag or the heart lung machine and allow blood to flow into the Fibrin Preparation Unit. Warning: Heavy blood flow from heart lung machine might cause bubbles.

After preparation, immediately remove the Fibrin Preparation Unit from the Processor Unit, turn it upside down and remove the fibrin syringe from the Fibrin Preparation Unit.

After preparation, store the fibrin syringe at room temperature. The fibrin syringe should be used within eight hours of preparation.

SPRAYPEN KIT / ENDOSCOPIC KIT - INSTRUCTIONS FOR USE

Vivostat Spraypen Kit / Endoscopic Kit contents

A sterile disposable single-use kit, which contains all the components needed to apply the sealant.

- **Spraypen Applicator/Endoscopic Applicator** - The application device is used to mix and apply the sealant.
- **pH Paper** - The pH paper is used to test adequate mixing of the fibrin.
- **pH 10** - The pH 10 solution activates fibrin polymerisation on spraying.
- **21 G Needle** - The 21 G needle is used to draw pH 10 buffer into the 1 ml syringe.
- **1 ml Syringe** - The 1 ml syringe is used to deliver pH 10 buffer via the application device. Caution: Avoid air bubbles when filling with pH10.
- **Batch Tracking labels (4 pcs.)** - Used for documentation in the patient journal.

Contraindications

Known hypersensitivity to batroxobin and tranexamic acid.

Cautions and Warnings

It is recommended that, whenever possible, Vivostat® sealant should be prepared from patient-derived blood or plasma.

Patient-derived sealant is for autologous use only – check patient details to ensure patient donor and recipient are the same! Discard any remaining sealant at the end of the procedure.

Note: When using blood or plasma derived from a source other than the patient, the risk of transmission of blood-borne infectious agents cannot be entirely excluded.

Spray is achieved using pressurised air. Caution should be used where there is a risk of air embolism. We recommend to use the spray mode "Jet No Air" in neurosurgery.

Direct topical application to the central nervous system of sufficiently high amounts of tranexamic acid (including application in fibrin sealants) may induce hyperexcitability, which may lead to convulsions. However, the concentration of tranexamic acid in Vivostat® (about 0.5mg/ml) has been shown to be between 6 and 12 times lower than that required to produce any noticeable effects in a rat study, where it was applied directly to the exposed spinal cord. Furthermore, clinical studies using Vivostat® in neurosurgery have identified no adverse effects related to the use of Vivostat®. However, the means by which hyperexcitability is induced by tranexamic acid and the importance of other factors when applied in fibrin sealants (e.g. individual variation, threshold in concentration/amount applied, formulation and release characteristics, mode and site of application) is poorly understood. Although indications are that the risk of inducing convulsions with Vivostat® is small, the possibility of inducing convulsions in susceptible individuals cannot be ruled out. Therefore, where direct contact of Vivostat® is made with the central nervous system, such as in spinal or neurosurgery, the operating staff should be aware of the risk of convulsions and be prepared to take appropriate corrective action if they should occur.

Not for intravascular use.

Dispose of properly after use.

Universal precautions for blood handling should always be used in the operation of this product. Dispose of all accessories that may contain biohazard materials, such as body fluids, using universal blood handling precautions.

Interactions

Drugs that reduce fibrinogen levels such as batroxobin and ancred may affect the preparation of the sealant. Administration of cefuroxime just prior to drawing blood may affect the coagulation abilities of the fibrin. The effect of acute thrombolytics such as streptokinase, urokinase and t-PA on the preparation and fibrinolytic stability of the sealant has not been evaluated and cannot be excluded.

Studies have been conducted to demonstrate that heparin, warfarin and aspirin therapy do not affect the preparation of fibrin using the Vivostat® system.

No other interactions are currently known.

Packaging, Storage and Shelf Life

Store between 4°C and 25°C.

Do not use after the expiry date on package labels.

Do not use if the packaging of any components appears open or damaged.

Product and Patient Tracking

A patient ID label should be placed on the fibrin syringe and not on the Fibrin Preparation Unit, as this will interfere with the process.

The label placed on the fibrin syringe should be checked prior to use for patient identification purposes.

Application Technique

- Dab the area of application until it is dry.
- In case of diffuse bleedings, lightly press a piece of gauze against the bleeding surfaces. Slowly remove the gauze and gradually apply the fibrin sealant to the application surface **in "High" mode** as soon as possible.
- The spray mode can be adjusted using the High/ Low button on the 300/350 AAU or the User Button on the 400/404 APL.
- In case of blockage of the Spraypen, cut off up to 8 mm of the tip of the Spraypen with a sharp scalpel.
- In case you change to a new Spraypen always remember to refill the pH10 syringe.

Cyclokapron® is a registered trade mark of Pharmacia

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Vivostat A/S