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VIVOSTAT® SYSTEM

Introduction

This PRF®-Preparation Kit / Co-Delivery Kit is part of a system and these instructions should be read in conjunction with the Vivostat® User Manuals and Quick Reference Cards.

The Vivostat® system comprises a Processor Unit, a disposable PRF® Preparation Kit for preparation of platelet rich fibrin, an Applicator Unit with a footswitch and a disposable Co-Delivery kit used to apply platelet rich fibrin and the desired substance (drug, stem cell etc.) during surgery.

Intended use

The Vivostat® system is a medical device used for the preparation and application of a platelet rich fibrin from whole blood. The Co-Delivery Kit has the ability to co-apply a fluid (substance) in a separate tube.

Choice of substance and its influence on the patient is the surgeon's responsibility.

Platelet rich fibrin has different scopes of application: Soft tissue regeneration, chronic wound healing, haemostasis, tissue sealing and tissue gluing.

PRF®-PREPARATION KIT – INSTRUCTIONS FOR USE

Vivostat PRF® Preparation Kit Contents

A sterile disposable single use kit, which contains all the components needed for the preparation of platelet rich fibrin.

- **PRF® Preparation Unit** - The PRF® Preparation Unit is a device that extracts platelets and fibrin from the blood into platelet rich fibrin. A Tubing Assembly is attached to the PRF® Preparation Unit used to fill the Preparation Unit.
- **Citrate/TA(PR®)** - a vial containing a sterile solution of citrate (anticoagulant) and 100 mg tranexamic acid (antifibrinolytic).
- **Codan Spike** - The dispensing spike is used to access the citrate vial and connected to the PRF® Preparation Unit.
- **Sugi swab** - The two non sterile Sugi swaps are used to clean the valve after the tubing assembly has been removed from the PRF® Preparation Unit.
- **pH4** - the pH4 syringe is loaded into the PRF® Preparation Unit prior to processing. After processing, this syringe contains the platelet rich fibrin solution and is referred to as the fibrin syringe.
- **Replacement Cap** - The Replacement Cap is used to protect the luer tip of the fibrin syringe.
- **Disinfection swab** - Used for disinfection of the rubber membrane on the Citrate/TA (PRF®)
- **Winged infusion set (18G)** - Used for blood donation connected to the tubing line on the Preparation Unit.
- **Batch Tracking labels (4 pcs.)** - Used for documentation in the patient journal.

Filling the PRF®-Preparation Unit with Blood

The Vivostat PRF® solution can only be prepared from patient-derived blood. We recommend to draw blood directly from the patient before any procedure is initiated.

Please notice! When the PRF® Preparation Unit contains the blood it must be handled carefully. Avoid shaking movements from side to side. Place the PRF® Preparation Unit in the Processor Unit and process it within one hour after drawing the blood.

After preparation, immediately remove the PRF® Preparation Unit from the Processor Unit, turn it upside down and remove the fibrin syringe from the PRF® Preparation Unit.

After preparation, store the fibrin syringe at room temperature. The fibrin syringe should be used within eight hours of preparation and be carefully inverted right before use.

CO-DELIVERY KIT - INSTRUCTIONS FOR USE

This Co-Delivery Kit is part of a system and these instructions should be read in conjunction with the Vivostat® User Manuals and Quick Reference Cards.

The Co-Delivery Kit only works together with the Applicator Units, APL 404.

The application of PRF® and and the additional substance can only be activated with the Vivostat footswitch as there is no pen button on the Co-Delivery Spraypen.

Vivostat® Co-Delivery Kit Contents

A disposable single-use kit, which contains all the components needed to apply the Platelet Rich Fibrin.

- **Co-Delivery Spraypen/Endoscopic Catheter** - Used to mix and apply the platelet rich fibrin with the additional substance.
- **pH Paper** - The pH paper is used to test adequate mixing of the Platelet Rich Fibrin.
- **pH 10** - The pH 10 solution activates the polymerisation of Platelet Rich Fibrin when spraying.
- **21 G Needle** - The 21 G needle is used to draw the pH10 Buffer into the 1 mL syringe.
- **1 mL Syringes** - The 1 ml syringe is used to deliver pH 10 buffer via the application device.
- **Co-Delivery Syringe** - For delivering a Platelet Rich Fibrin compatible substance via the Co-Delivery Spraypen.
Caution: Avoid air bubbles in both syringes when filling with pH-10 and the Platelet Rich Fibrin compatible substance.
N.B. You can not fill more substance in the Co-Delivery syringe than is found in the fibrin syringe. See illustration on the last page.
- **Replacement Cap** - The Replacement Cap is placed on the tip of the Co-Delivery syringe if it is not used immediately after filling.

- **Batch Tracking labels (4 pcs.)** - Used for documentation in the patient journal.

Contraindications

Known hypersensitivity to batroxobin and tranexamic acid.

Known hypersensitivity to the chosen additional substance.

Cautions and Warnings

It is recommended to prepare Vivostat® platelet rich fibrin from patient-derived whole blood only.

Patient-derived platelet rich fibrin is for autologous use only – check patient details to ensure patient donor and recipient is the same!

Spray is achieved using pressurised air. Caution should be used where there is a risk of air embolism.

Direct topical application to the central nervous system of sufficiently high amounts of tranexamic acid (including application in fibrin sealants) may induce hyperexcitability, which may lead to convulsions. However, the concentration of tranexamic acid in Vivostat® (about 0.5mg/mL) has been shown to be between 6 and 12 times lower than that required to produce any noticeable effects in a rat study, where it was applied directly to the exposed spinal cord. Furthermore, clinical studies using Vivostat® in neurosurgery have identified no adverse effects related to the use of Vivostat®. However, the means by which hyperexcitability is induced by tranexamic acid and the importance of other factors when applied in fibrin sealants (e.g. individual variation, threshold in concentration/amount applied, formulation and release characteristics, mode and site of application) are poorly understood. Therefore, although indications are that the risk of inducing convulsions with Vivostat® is small, the possibility of inducing convulsions in susceptible individuals cannot be ruled out. Therefore, where direct contact of the Vivostat® solution is made with the central nervous system, such as in spinal or neurosurgery, the operating staff should be aware of the possibility of convulsions and be prepared to take appropriate corrective action if they should occur.

Choice of additional substance and its influence on the platelet rich fibrin is the surgeon's responsibility.

Not for intravascular use.

Dispose properly after use.

Universal precautions for blood handling should always be used in the operation of this product. Dispose of all accessories that may contain biohazardous materials, such as body fluids, using universal blood handling precautions.

Interactions

Drugs that reduce fibrinogen levels such as batroxobin and ancrod may affect the preparation of the platelet rich fibrin. **Administration of cefuroxime just prior to drawing blood may affect the coagulation abilities of fibrin.** The effect of acute thrombolytics such as streptokinase, urokinase and t-PA on the preparation and fibrinolytic stability of the

platelet rich fibrin has not been evaluated and cannot be excluded.

Studies have been conducted to demonstrate that heparin, warfarin and aspirin therapy do not affect the preparation of fibrin using the Vivostat® System.

Platelet aggregation and release of platelet-derived growth factors may be affected for up to 7 days after ingestion of analgesics containing acetylsalicylic acid. Other anti-platelets drugs (analgesics) may also have the same effect. The consequence of this for performance of Vivostat PRF® is not known currently.

Trombocytopeni in the patient may cause lower contents of platelets than expected.

No other interactions are currently known.

Packaging, Storage and Shelf Life

Store between 4°C and 25° C.

Do not use after the expiry date on package labels

Do not use if the packaging of any components appears to be open or damaged.

Product and Patient Tracking

Labels should be placed on the fibrin syringe and not on the PRF®-Preparation Unit, as this will interfere with the process.

Labels placed on the fibrin syringe should be checked prior to use for patient identification purposes.

Application Technique

- Dab the area of application until it is dry.
- The spray mode on the APL 404 Applicator can be adjusted on the User Button. We recommend to use the spraymode „low“ for Co-Delivery applications.
- In case of blockage of the Spraypen, cut off upto 8mm of the tip of the Spraypen with a sharp scalpel.
- If you change to a new Spraypen always remember to refill the pH-10 syringe.

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Vivostat A/S