PRF®-PREPARATION KIT & SPRAYPEN KIT

(SPRAYPEN KIT & ENDOSCOPIC KIT)

Manufactured by:

Vivostat A/S
Borupvang 2
3450 Allerød
Denmark
www.vivostat.com

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THE VIVOSTAT® SYSTEM

Introduction

This PRF® Preparation Kit and Spraypen Kit / Endoscopic Kit is part of a system and these instructions should be read in conjunction with the Vivostat® User Manuals and Quick Reference Cards.

The Vivostat® System comprises: A Processor Unit, a disposable PRF® Preparation Kit for preparation of platelet rich fibrin, an Applicator Unit and a disposable Application Kit used to apply platelet rich fibrin during surgery.

Endoscopic Applicator Handle: The Endoscopic Applicator Handle is a reusable device, used to insert and manipulate the Endoscopic catheter.

Intended use

The Vivostat® system is a medical device used for the preparation and application of a platelet rich fibrin from whole blood.

Platelet rich fibrin has different scopes of application: Soft tissue regeneration, chronic wound healing, haemostasis, tissue sealing and tissue gluing.

PRF® PREPARATION KIT – INSTRUCTION FOR USE

Vivostat PRF® Preparation Kit contents

A sterile disposable single use kit, which contains all the components needed for the preparation of platelet rich fibrin.

- **PRF® Preparation Unit** – The PRF® Preparation Unit is a device that extracts platelets and fibrin from the blood into platelet rich fibrin. A tubing assembly is attached to the PRF® Preparation Unit.
- **Citrate/TA (PRF®)** – a vial containing a sterile solution of citrate (anticoagulant) and 100 mg tranexamic acid (antifibrinolytic).
- **Spike** – The dispensing spike is used to access the Citrate Vial and connect to the PRF® Preparation Unit.
- **Sugi swap** – The two non sterile Sugi swaps are used to clean the valve after the tubing assembly has been removed from the PRF® Preparation Unit.
- **pH4** – the pH4 syringe is loaded into the PRF® Preparation Unit prior to processing. **After processing, this syringe contains the platelet rich fibrin solution and is referred to as the fibrin syringe.**
- **Replacement Cap** – The Replacement Cap is used to protect the luer tip of the fibrin syringe.
- **Disinfection swap** – Used for disinfection of the rubber membrane on the Citrate/TA (PRF®).
- **Winged infusion set (18G)** – Used for blood donation connected to the tubing line at the Preparation Unit.
- **Removable Batch Tracking labels (4 pcs.)** – Used for documentation in the patient journal.

Filling the PRF® Preparation Unit with Blood

Vivostat® platelet rich fibrin can only be prepared from patient derived donation of blood. We recommend to draw blood directly from the patient before any procedure is initiated.

Please notice! When the PRF® Preparation Unit contains the blood it must be handled carefully. Avoid shaking movements from side to side. The PRF® Preparation Unit containing the blood must be processed within one hour after collecting the blood.

After preparation, immediately remove the PRF® Preparation Unit from the Processor Unit, turn it upside down and remove the fibrin syringe from the PRF® Preparation Unit.

After preparation, store the fibrin syringe at room temperature. The fibrin syringe should be used within eight hours of preparation and be carefully inverted right before use.

SPRAYPEN KIT / ENDOSCOPIC KIT – INSTRUCTIONS FOR USE

Vivostat® Spraypen Kit / Endoscopic Kit contents

A sterile disposable single use kit, which contains all the components needed to apply the platelet rich fibrin.

- **Spraypen Applicator/Endoscopic Applicator** – The application device is used to mix and apply the platelet rich fibrin.
- **pH Paper** – The pH paper is used to test adequate mixing of the platelet rich fibrin.
- **pH 10** – The pH 10 solution activates the polymerisation of the platelet rich fibrin on spraying.
- **21 G Needle** – The 21 G needle is used to draw pH 10 buffer into the 1 ml syringe.
- **1 ml Syringe** – The 1 ml syringe is used to deliver pH 10 buffer via the application device. Caution: Avoid air bubbles when filling with pH10.
- **Batch Tracking labels (4 pcs.)** – Used for documentation in the patient journal.

Contraindications

Known hypersensitivity to batroxobin and tranexamic acid.

Cautions and Warnings

It is recommended to prepare Vivostat® platelet rich fibrin from patient derived whole blood only.

Patient derived platelet rich fibrin is for autologous use only – check patient details to ensure patient donor and recipient are the same!

Spray is achieved using pressurised air. Caution should be used where there is a risk of air embolism. We recommend to use the spray mode "Jet No Air" in neurosurgery.

Direct topical application to the central nervous system of sufficiently high amounts of tranexamic acid

(including application in fibrin sealants) may induce hyperexcitability, which may lead to convulsions. However, the concentration of tranexamic acid in Vivostat® (about 0.5mg/ml) has been shown to be between 6 and 12 times lower than that required to produce any noticeable effects in a rat study, where it was applied directly to the exposed spinal cord. Furthermore, clinical studies using Vivostat® in neurosurgery have identified no adverse effects related to the use of Vivostat®. However, the means by which hyperexcitability is induced by tranexamic acid and the importance of other factors when applied in fibrin sealants (e.g. individual variation, threshold in concentration/amount applied, formulation and release characteristics, mode and site of application) is poorly understood. Although indications are that the risk of inducing convulsions with Vivostat® is small, the possibility of inducing convulsions in susceptible individuals cannot be ruled out. Therefore, where direct contact of Vivostat® is made with the central nervous system, such as in spinal or neurosurgery, the operating staff should be aware of the risk of convulsions and be prepared to take appropriate corrective action if they should occur.

Not for intravascular use.

Dispose of properly after use.

Universal precautions for blood handling should always be used in the operation of this product. Dispose of all accessories that may contain biohazard materials, such as body fluids, using universal blood handling precautions.

Interactions

Drugs that reduce fibrinogen levels such as batroxobin and ancred may affect the preparation of the platelet rich fibrin. Administration of cefuroxime just prior to drawing blood may affect the coagulation abilities of fibrin. The effect of acute thrombolytics such as streptokinase, urokinase and t-PA on the preparation and fibrinolytic stability of the platelet rich fibrin has not been evaluated and can not be excluded.

Studies have been conducted to demonstrate that heparin and warfarin therapy do not affect the preparation of fibrin using the Vivostat® System.

Platelet aggregation and release of platelet-derived growth factors may be affected for up to 7 days after ingestion of analgesics containing acetylsalicylic acid. Other anti-platelets drugs (analgesics) may also have the same effect. The consequence of this for performance of Vivostat PRF® is not known currently.

Trombocytopeni in the patient may cause lower contents of platelets than expected.

No other interactions are currently known.

Packaging, Storage and Shelf Life

Store between 4°C and 25°C.

Do not use after the expiry date on package labels.

Do not use if the packaging of any components appears open or damaged.

Product and Patient Tracking

A patient ID label should be placed on the fibrin syringe and not on the PRF®-Preparation Unit, as this will interfere with the process.

The label placed on the fibrin syringe should be checked prior to use for patient identification purposes.

Application Technique

- Dab the area of application until it is dry.
- In case of diffuse bleedings, lightly press a piece of gauze against the bleeding surfaces. Slowly remove the gauze and gradually apply the fibrin sealant to the application surface **in "High" mode** as soon as possible.
- The spray mode can be adjusted using the High/Low button on the 300/350 AAU or the User Button on the 400/404 APL.
- In case of blockage of the Spraypen, cut off up to 8 mm of the tip of the Spraypen with a sharp scalpel.
- In case you change to a new Spraypen always remember to refill the pH10 syringe

Cyclokapron® is a registered trade mark of Pharmacia

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