Vivostat® Split Kit

Instructions for use

Manufactured by:

Vivostat A/S
Borupvang 2
3450 Alleroed
Denmark
www.vivostat.com

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SPLIT KIT - INSTRUCTION FOR USE

Introduction

This Split Kit is part of a system and these instructions should be read in conjunction with the Vivostat® system User Manual and Instruction For Use (IFU) found in the kit used to prepare PRF®/Fibrin.

Intended Use

Autologous Platelet Rich Fibrin (PRF®) can be used for healing chronic wounds. The intended use of the Split Kit is to freeze the PRF®/Fibrin solution with the intention to treat the patient more than one time from one blood donation.

Vivostat Split Kit Components

- Syringe with transfer tube attached
- One Cap - to close the tip of the syringe after transferring PRF® or Fibrin the syringe
- Piston Rod: Used to transfer PRF®/Fibrin from one syringe to another

SYMBOLS

The Vivostat® system is CE marked in accordance with the provision of the EC Medical Device Directive 93/42/EEC

Note : See Instructions for Use / User Manual

Single use (Do not reuse)

REF

Catalog / Order number / Item No.

Lot/Batch number

Expiration date (Use By)

Sterilised using ethylene oxide (EtO)

Contraindications

Known hypersensitivity to batroxobin or Tranexamic Acid

Cautions and Warnings

Patient derived fibrin and platelet rich fibrin is for autologous use only – check patient details to ensure patient donor and recipient are the same! Do not re-freeze the PRF®/Fibrin solution!

It is important for the efficiency of the PRF®/Fibrin product that time of thawing is observed. It is also particularly important to use the PRF® solution within 20 minutes after thawing. The coagulation of PRF®/Fibrin may in some cases be prolonged.

For practical reasons it is recommended that thawing of PRF® takes place near the operat-

ing/treatment-room.

Packaging, Storage and Shelf Life

Do not use after expiry date on single use packaging or package labels.

Do not use if the packaging of any components appears open or damaged.

The frozen product must be used within 3 weeks from the time it was frozen.

Product and Patient Tracking

Labels with patient identification and date of freezing should as a minimum be placed on the syringe before freezing.

Before thawing, labels, must be checked prior to use for correct patient identification purposes.

Filling the Syringe

1. Open the Split Kit containing a syringe, one cap and a piston rod.
2. Apply a patient identification label to the empty syringe.
3. Remove all air above the PRF[®]/Fibrin solution in the filled syringe by pushing the plunger upwards in vertical direction with the piston rod.
4. Attach the filled syringe to the empty syringe by means of the preattached transfer tube.
5. Use the piston rod to push the desired volume of PRF[®]/Fibrin from the filled syringe to the empty syringe
6. Remove the transfer tube from the syringe (you may use a pean or similar) by turning and pulling upwards.
7. Attach the cap to the syringe. Correct placement will be confirmed by a click.
8. Place the syringe with patient identification in a freezer with temperature control, between -18°C and -25°C.

Thawing of PRF[®]/Fibrin

1. Thaw the PRF[®]/Fibrin syringe at room temperature for 40 minutes.
2. Mix the contents of the syringe carefully by turning it upside down about 10 times.
3. Spray with PRF[®]/Fibrin within 20 minutes from thawing. The time of coagulation may be some what slower than that of freshly prepared PRF[®]/Fibrin.
4. Thawed PRF[®]/Fibrin may not appear as firm as PRF[®]/Fibrin that has not been frozen.
5. Discard any remaining PRF[®]/Fibrin after use. **Do not re-freeze!**

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